




**MICRO
OPHTHALMIC
KNIFE**



LANCE TIP/ SIDE PORT BLADE

Color: Yellow

FOR SCLERAL INCISION / SIDE PORT

Product Code	Size		Type
OPLTB1	15	24G	Straight / Angled
OPLTB2	30	24G	Straight / Angled
OPLTB3	45	24G	Straight / Angled

For ordering Angled Knife after code add "AN"



CRESCENT BLADE

Color: Pink

FOR MAKING TUNNEL

Product Code	Size / Description		
OPCB1	Angled 45	Bevel up	20G
OPCB2	Angled 45	Bevel down	20G
OPCB3	Angled 45	Bevel up	19G



KERATOME (SLIT) BLADE

Color: Blue

For Entry into Anterior Chamber

Angled 45 Bevel up (Pointed Tip) 19 G

Also known as clear corneal blade

Double Bevel

Size	Code
2.8mm	KBDB28
3.8mm	KBDB38

Product Code	Code
1.8mm	OPKB18
2.2mm	OPKB22
2.5mm	OPKB25
2.6mm	OPKB26
2.8mm	OPKB28
3.0mm	OPKB30
3.2mm	OPKB32
3.5mm	OPKB35




EXTENSION KERATOME BLADE

Color: Green

FOR IOL ENLARGING

Blunt tip 45 Angled Bevel up

Product Code	Size	
OPKE1	5.1mm	16G
OPKE2	5.2mm	16G
OPKE3	5.5mm	15G
OPKE4	6.0mm	15G



MVR BLADE

Color: Pink

FOR WATER TIGHT SEALING INCISION

Product Code	Size / Description	
OPMV19	19G	Straight/ Angled
OPMV20	20G	Straight/ Angled
OPMV24	24G	Straight/ Angled

For ordering Angled Knife after code add "AN"