

## Green Max “Regular” Fiber Optic Reusable Blades

Product Model : **LSC22**



**Category: Medical Equipment, Laryngoscope Set & Spares,**

### Summary :

Utilizing advanced LED light technology, **GWS Green Max Fiber Optic laryngoscopes** deliver superior illumination for safer intubations. Reusable Fiber Optic blades are manufactured with 100% high quality stainless steel to ensure durability and long life. Each blade is MATT finishes to reduce glare and enhance ease of intubation for optimal patient safety.

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Enlarged 3.6mm diameter fiber optic bundle.
- Engineered with a removable fiber optic bundle for easy replacement of damaged bundles.
- Autoclavable Light Guide, tested to withstand over 4000 cycles.
- Matt Finish to avoid any reflection during intubation.
- Compatible to all ISO Handles 7376-3
- Available in Macintosh & Miller Blades.

Size	Item Code
MAC0	LSC22.1
MAC1	LSC22.2
MAC2	LSC22.3
MAC3	LSC22.4
MAC4	LSC22.5
MAC5	LSC22.6
MILL00	LSC22.7
MILL0	LSC22.8
MILL1	LSC22.9
MILL2	LSC22.10
MILL3	LSC22.11
MILL4	LSC22.12