

Kidney Tray, Stainless Steel

Product Model : **HH-7**



Category: Stainless Steel Hospital Hollowware, General Medical Devices,

Summary :

Specification :

Kidney Trays (Stainless Steel) without cover					
Approx. Size		Quality / Code			
In (mm)	In (Inches)	Export Quality/ Thickness: 0.35 mm	Deluxe Quality/ Thickness: 0.5 mm	Superior Quality/ Grade 304, 24G / Thickness: 0.5 mm	Super Deluxe Quality/ Grade 304, 24G /Thickness: 0.63 mm

150 mm	6"	HH-7.1.6	HH-7.2.6	HH-7.3.6	HH-7.4.6
200 mm	8"	HH-7.1.8	HH-7.2.8	HH-7.3.8	HH-7.4.8
250 mm	10"	HH-7.1.10	HH-7.2.10	HH-7.3.10	HH-7.4.10
300 mm	12"	HH-7.1.12	HH-7.2.12	HH-7.3.12	HH-7.4.12
250 mm - Shallow type (25x15x4cm)	10"	---	---	---	HH-7.4.S10

Note: Kidney Tray with Cover also available. Add Code "C" before the above codes.

Legacy of Excellence in Health Care..



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com