

## Kidney Tray, Stainless Steel

Product Model : **HH-7**



**Category: Stainless Steel Hospital Hollowware, Medical Equipment,**

**Summary :** Kidney Tray, Stainless Steel is a durable and hygienic medical tool designed for use in hospitals, clinics, and medical facilities. Made from high-grade stainless steel, this tray is corrosion-resistant, easy to clean, and maintains its shine even after repeated use. Its ergonomic shape and smooth edges ensure patient comfort and ease of handling for medical professionals. Ideal for holding surgical instruments, specimens, and other medical items, this kidney tray is an essential component of any medical environment.

**Specification :**

Kidney Trays (Stainless Steel) without cover	
Approx. Size	Quality / Code

In (mm)	In (Inches)	Export Quality/ Thickness: 0.35 mm	Deluxe Quality/ Thickness: 0.5 mm	Superior Quality/ Grade 304, 24G / Thickness: 0.5 mm	Super Deluxe Quality/ Grade 304, 24G  /Thickness: 0.63 mm
150 mm	6"	HH-7.1.6	HH-7.2.6	HH-7.3.6	HH-7.4.6
200 mm	8"	HH-7.1.8	HH-7.2.8	HH-7.3.8	HH-7.4.8
250 mm	10"	HH-7.1.10	HH-7.2.10	HH-7.3.10	HH-7.4.10
300 mm	12"	HH-7.1.12	HH-7.2.12	HH-7.3.12	HH-7.4.12
250 mm - Shallow type (25x15x4cm)	10"	---	---	---	HH-7.4.S10
<i>Note: Kidney Tray with Cover also available. Add Code "C" before the above codes.</i>					

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com