

## Kidney Tray, Stainless Steel

Product Model : **HH-7**



**Category: Stainless Steel Hospital Hollowware, Medical Equipment,**

**Summary :** Kidney Tray, Stainless Steel is a durable and hygienic medical tool designed for use in hospitals, clinics, and medical facilities. Made from high-grade stainless steel, this tray is corrosion-resistant, easy to clean, and maintains its shine even after repeated use. Its ergonomic shape and smooth edges ensure patient comfort and ease of handling for medical professionals. Ideal for holding surgical instruments, specimens, and other medical items, this kidney tray is an essential component of any medical environment.

**Legacy of Excellence in Healthcare...**

✉ info@gwsmed.com, gwssurgicals@gmail.com

☎ +91 124 4915548

🌐 www.gwsmed.com



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

Kidney Trays (Stainless Steel) without cover					
Approx. Size		Quality / Code			
In (mm)	In (Inches)	Export Quality/ Thickness: 0.35 mm	Deluxe Quality/ Thickness: 0.5 mm	Superior Quality/ Grade 304, 24G / Thickness: 0.5 mm	Super Deluxe Quality/ Grade 304, 24G  /Thickness: 0.63 mm
150 mm	6"	HH-7.1.6	HH-7.2.6	HH-7.3.6	HH-7.4.6
200 mm	8"	HH-7.1.8	HH-7.2.8	HH-7.3.8	HH-7.4.8
250 mm	10"	HH-7.1.10	HH-7.2.10	HH-7.3.10	HH-7.4.10
300 mm	12"	HH-7.1.12	HH-7.2.12	HH-7.3.12	HH-7.4.12
250 mm - Shallow type (25x15x4cm)	10"	---	---	---	HH-7.4.S10
Note: Kidney Tray with Cover also available. Add Code "C" before the above codes.					