


5+3 Twin

Product Model : **OTE53**



Category: Medical Equipment, Operation Light, Elite Series Operation Light,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

Technical Data	
Intensity	1,60,000 + 1,00,000 Lux \pm 10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Reduction Index	93 RA
LED Life	?50000 Hrs.
Number of LED	105 + 63
Diameter of Light	600mm x 2
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional