


## 5+5 Twin

Product Model : **OTE55**



**Category: Medical Equipment, Operation Light, Elite Series Operation Light,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

### 5+5 OT Light

#### Technical Data

Intensity	1,60,000 + 1,60,000 Lux $\pm$ 10%
Size of Light Field	12-30 cm
Colour Temperature	4000-5000 K
Colour Reduction Index	93 RA
LED Life	?50000 Hrs.
Number of LED	105 + 105
Diameter of Light	600mm x 2
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional