

Electric OT Table Advanced

Product Model : **OT11**



Category: General Medical Devices, Operation Table,

Specification :

? SS-304 | ? Surgeons Panel | ? Electric (remote operated) | ? C Arm Compatible | ? Floor Locking | ? Zero Positioning

- Fully electrically operated table, five sections, radio-translucent table with eccentric pillars and Die compressed base.
- Interchangeable head and leg sections (manually)
- All positions operated through cable remote.
- Electric floor locking. (Optional)
- Inbuilt kidney bridge.
- Auto leveling (zero function)
- We are providing two handsets. The additional fixed control panel provided on the tabletop, which can be used in case of failure of the main handset.
- Indication of battery back provided in the handset.

• The handset provides dual memory functions which can be set and reset.

Top dimension	L 1900mm x W 522mm
Height adjustment	750mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney elevator	130mm
Backrest (up/down)	60° / 25°
Leg rest (up/down)	15° / 90°
Headrest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains Voltage	AC 220V ±10% 50Hz
Patient weight capacity	200 kgs.

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicalsgmail.com



+91-124-4915548



www.gwsmed.com