


## Electric OT Table Advanced

Product Model : **OT11**



**Category: Medical Equipment, Operation Table,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

? SS-304 | ? Surgeons Panel | ? Electric (remote operated) | ? C Arm Compatible | ? Floor Locking | ? Zero Positioning

- Fully electrically operated table, five sections, radio-translucent table with eccentric pillars and Die compressed base.
- Interchangeable head and leg sections (manually)
- All positions operated through cable remote.
- Electric floor locking. (Optional)
- Inbuilt kidney bridge.
- Auto leveling (zero function)
- We are providing two handsets. The additional fixed control panel provided on the tabletop, which can be used in case of failure of the main handset.
- Indication of battery back provided in the handset.
- The handset provides dual memory functions which can be set and reset.

Top dimension	L 1900mm x W 522mm
Height adjustment	750mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney elevator	130mm
Backrest (up/down)	60° / 25°
Leg rest (up/down)	15° / 90°
Headrest (up/down)	20° / 60°
Power Supply	24V DC
Battery Backup	2-3 hours (optional)
Mains Voltage	AC 220V ±10% 50Hz
Patient weight capacity	200 kgs.