

## Electro Hydraulic OT Table

Product Model : **OT13**



**Category: Medical Equipment, Operation Table,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

### Features:

- Tabletop is made with Stainless Steel grade 304 to avoid rust with a canceled pump.
- The five-section stainless steel tabletop is equipped with removable, antistatic mattresses.
- All the positions are realized by a mechanical driver.
- The table is furnished with a hard, robust stainless steel shell which provides the sense of dirt-free & sanitized operations.
- Precision-made hydraulic system with single lever action hydraulic life controlled by a foot pedal.
- Five sections stainless steel top gives anatomic adaptation to the patient body.
- Inbuilt kidney bridge and large perineal cut out for drainage.
- Manual Floor locking system.

### Technical Data:

Top dimension	L 1900mm x W 520mm
Height adjustment	762mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney elevator	150mm
Back rest (up/down)	80° / 25°
Leg rest (up/down)	15° / 90°
Head rest (up/down)	20° / 60°
Patient weight capacity	200 kgs.