

## **Electro Hydraulic OT Table**

Product Model : **OT13**



**Category: General Medical Devices, Operation Table,**

**Specification :**

### **Features:**

- Tabletop is made with Stainless Steel grade 304 to avoid rust with a canceled pump.
- The five-section stainless steel tabletop is equipped with removable, antistatic mattresses.
- All the positions are realized by a mechanical driver.
- The table is furnished with a hard, robust stainless steel shell which provides the sense of dirt-free & sanitized operations.
- Precision-made hydraulic system with single lever action hydraulic life controlled by a foot pedal.
- Five sections stainless steel top gives anatomic adaptation to the patient body.
- Inbuilt kidney bridge and large perineal cut out for drainage.
- Manual Floor locking system.

**Technical Data:**

Top dimension	L 1900mm x W 520mm
Height adjustment	762mm - 1000 mm
Trendelenburg / Reverse	30° / 25°
Lateral Tilt	20° / 20°
Kidney elevator	150mm
Back rest (up/down)	80° / 25°
Leg rest (up/down)	15° / 90°
Head rest (up/down)	20° / 60°
Patient weight capacity	200 kgs.

*Legacy of Excellence in Health Care..*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com