

Endotracheal Tube Plain & Cuffed

Product Model : **MD12**



Category: Medical Equipment, Medical Disposables, Anaesthesia,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Designed to be inserted into the trachea through the oral/nasal cavity to secure the airway and maintain the gas exchange in the lungs.
- Made from non-toxic, non-irritant medical grade PVC
- Radio opaque line provided throughout the length for X-ray visualization
- Universal 15mm Connector at Proximal End
- The tip is suitable for both Nasal & Oral intubation
- Cuffed Sizes (mm) : 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10mm

Endotracheal Tube Cuffed - Sizes/ Codes

Endotracheal Tube Cuffed - Sizes/ Codes								
Size	2.0mm	2.5mm	3.0mm	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm
Code	MD12.1	MD12.2	MD12.3	MD12.4	MD12.5	MD12.6	MD12.7	MD12.8
Size	6.5mm	7.0mm	7.5mm	8.0mm	8.5mm	9.0mm	9.5mm	10mm
Code	MD12.9	MD12.10	MD12.11	MD12.12	MD12.13	MD12.14	MD12.15	MD12.16

Packing

- ETO Sterile, individually packed in a peelable pouch pack.
- Master Carton of 200 pcs. (Apprx.)

Gross Weight:

6.7 kg (300 pcs.)

Master Carton:

660 x 320 x 380mm (Size:3.0 to 4.5mm)

Gross Weight:

12.5 kg (300 pcs.)

Master Carton:

660 x 380 x 50mm (Size:5.0 to 9.5mm)