

Endotracheal Tube Plain & Cuffed

Product Model : **MD12**



Category: Medical Equipment, Medical Disposables, Anaesthesia,

Specification :

- Designed to be inserted into the trachea through the oral/nasal cavity to secure the airway and maintain the gas exchange in the lungs.
- Made from non-toxic, non-irritant medical grade PVC
- Radio opaque line provided throughout the length for X-ray visualization
- Universal 15mm Connector at Proximal End
- The tip is suitable for both Nasal & Oral intubation
- Cuffed Sizes (mm) : 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10mm

Endotracheal Tube Cuffed - Sizes/ Codes

Size	2.0mm	2.5mm	3.0mm	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm
Code	MD12.1	MD12.2	MD12.3	MD12.4	MD12.5	MD12.6	MD12.7	MD12.8