

Examination Table with Height Adjustment, Hydraulic

Product Model : **HF145**



Category: hospital furniture, General Medical Devices, Examination Couch & Tables,

Specification :

- Two Section Top
- MS Framework, pre-treated & epoxy powder coated.
- Adjustable Height by Hydraulic Pump
- Weight Load Capacity 180-200Kgs
- Height Adjustment(Lo-Hi):- 51cms to 82cms
- Product Dimensions:- 72? x 24? x 20?
- Also Available in Electric:- **EHF145**

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwsurgicals@gmail.com  +91-124-4915548  www.gwsmed.com