

## Sterilization Tray with Silicone Mat, Stainless Steel

Product Model : **OST11**



**Category: Stainless Steel Hospital Holloware, Medical Equipment, Sterilization Trays,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Sterilization Tray, used for ophthalmic instruments sterilizing.
- Made with Stainless Steel Material.
- Silicone Mat inside.
- Press button lock for opening and closing

Item Code	Type	Size(cms) LxWxH	Net Weight
OST11.01	Small	20 x 6.5 x 2.5	230 Grams
OST11.02	Medium	21 x 12.5 x 2.5	400 Grams
OST11.03	Large	on request	on request