

## Sterilization Tray with Silicone Mat, Stainless Steel

Product Model : **OST11**



**Category: Stainless Steel Hospital Hollowware, General Medical Devices, Sterilization Trays,**

### Specification :

- Sterilization Tray, used for ophthalmic instruments sterilizing.
- Made with Stainless Steel Material.
- Silicone Mat inside.
- Press button lock for opening and closing

| Item Code | Type   | Size(cms) LxWxH | Net Weight |
|-----------|--------|-----------------|------------|
| OST11,01  | Small  | 20 x 6.5 x 2.5  | 230 Grams  |
| OST11,02  | Medium | 21 x 12.5 x 2.5 | 400 Grams  |
| OST11,03  | Large  | on request      | on request |

*Legacy of Excellence in Health Care...*

Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

info@gwsmed.com, gwssurgicals@gmail.com +91-124-4915548 www.gwsmed.com