

## Silicone Artificial Resuscitator-Adult (with POP off Valve)

Product Model : **AN101**



**Category: Anaesthesia Products & Equipment, Medical Equipment, Silicone Artificial Resuscitators (Ambu Bag),**

### Specification :

- Capacity: 1600 ml.
- Oxygen reservoir bag: 2600 ml.
- Resuscitation set with 60cm H<sub>2</sub>O pop of valve.
- Disposable O<sub>2</sub> tubing 2 mtr.
- Reusable & latex free.
- Silicone mask: 4 no.
- Color: Blue/Green/White
- Packing: Zip pouch packing.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com