

## Silicone Artificial Resuscitator-Adult (with POP off Valve)

Product Model : **AN101**



**Category: Anaesthesia Products & Equipment, Medical Equipment, Silicone Artificial Resuscitators (Ambu Bag),**

**Legacy of Excellence in Healthcare...**

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Capacity: 1600 ml.
- Oxygen reservoir bag: 2600 ml.
- Resuscitation set with 60cm H2O pop of valve.
- Disposable O2 tubing 2 mtr.
- Reusable & latex free.
- Silicone mask: 4 no.
- Color: Blue/Green/White
- Packing: Zip pouch packing.