

## Black Resuscitator Non Autoclavable (Adult)

Product Model : **AN-108**



**Category: General Medical Devices, Anaesthesia Products & Equipment, Artificial Resuscitators(Ambu Bag) Black Rubber,**

### Specification :

- Double Layer Non-Autoclavable Bellow.
- Black Rubber Super Face Mask Size: 4 no.
- Disposable & Single Use
- Zip Pouch Packing
- Color-Black
- Packing: Zip Pouch Packing

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com