


## Five Functional ICU Bed Electric (ABS Panel & ABS Side Railing)

Product Model : **HF-101**



**Category:** hospital furniture, Medical Equipment, ICU Beds,

**Specification :**

- Advance Linear Actuator System [01pc. Motor + 01pc. Control Box + 01pc. Nurse Control System (ACP)]
- Four Individual Side Railing, ABS Moulded.
- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 460-700H mm
- Functions by Hand Remote: (a) Hi-low (b)Backrest elevation (c)Kneerest elevation (d)Trendlenburg / Reverse Trendlenburg
- Functions by Nurse Control System (ACP): (a)Hi-low with locking function for handset. (b)Backrest elevation with locking function for handset. (c)Kneerest elevation with locking function for handset. (d)Trendlenburg / Reverse Trendlenburg with locking function. (e)Electric CRP for Zero Positioning of Bed.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com