

RESPO-Infant Resuscitator

Product Model : **RESP01**



Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,

Specification :

- Used on delivery room, baby wards, NICU that the baby's maximum weight up to 10KG
- Inside an alarm system of differential pressure. When alarming, the safety valve starts automatically and produces the mixture gas continuously to ensure the safety of patients
- The range of flow adjustment is 5LPM~15LPM.
- The MAX-P and PIP are safe, controlled and accurate, PEEP is constant and accurate.

Technical Specifications:-

- Pressure gauge: -2~10kPa (-20~100cmH₂O).
- MAX-P: @0~15LPM 76kPa (60cmH₂O).
- PIP: @5~15LPM 0.2~5.7kPa (2~57cmH₂O).

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com