

RESPO(PRO)-Infant Resuscitator

Product Model : **RESP02**



Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Used in delivery room, baby wards, NICU that the baby's maximum weight up to 10KG
- Inside an alarm system of differential pressure. When alarming, the safety valve starts automatically and produces the mixture gas continuously to ensure the safety of patients
- Adjust the Oxygen concentration from 21% to 100% optionally
- The range of flow adjustment is 5LPM~15LPM.
- The MAX-P and PIP are safe, controlled and accurate, the PEEP is constant and accurate.

Technical Specifications:-

- Pressure gauge: -2~10kPa (-20~100cmH₂O).
- MAX-P: @0~15LPM 6kPa (60cmH₂O).
- PIP: @5~15LPM 0.2~5.7kPa (2~57cmH₂O).
- PEEP: @5~15LPM 0.03~2.3kPa (0.3~23cmH₂O).