

## **RESPO-CARE | Infant Bubble CPAP System**

Product Model : **RESP05**



**Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,**

### **Specification :**

GWS-RESP05 Bubble CPAP unit is full featured NICU equipment to meet neonatal critical care requirements for New-born baby those with spontaneous breathing infants requiring respiratory support. Complete system provides perfect delivery of FIO<sub>2</sub> with inbuilt Air / Oxygen Blender, Bubble generator with pressure manifold, heated breathing circuit and inbuilt medical grade Air compressor.

- **Servo Humidifier** Compatible for neonate to adult patients with Invasive & Non-invasive modes.
- **Blender** Compact and unique Air / Oxygen blender to maintain perfect FIO<sub>2</sub> from 21% to 100%. Oxygen flow directly can be connected to central pipe line or through flow meter regulator easy to use
- **Air Compressor** Low Noise Medical grade oil free low pressure Air compressor with bacteria filter, Maximum flow 15L/minute, High pressure relief value protection with constant air flow to the blender
- **Bubble Generator** CPAP Bubble generator with adjustable CPAP from 3 to 10 cmH<sub>2</sub>O, water chamber ease to fill and remove mounted on side of the stand

\*Accessories shown in the image may not be part of the standard machine.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com