

RESPO-CARE | Infant Bubble CPAP System

Product Model : **RESP05**



Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

GWS-RESP05 Bubble CPAP unit is full featured NICU equipment to meet neonatal critical care requirements for New-born baby those with spontaneous breathing infants requiring respiratory support.

Complete system provides perfect delivery of FIO₂ with inbuilt Air / Oxygen Blender, Bubble generator with pressure manifold, heated breathing circuit and inbuilt medical grade Air compressor.

- **Servo Humidifier**

Compatible for neonate to adult patients with Invasive & Non-invasive modes.

- **Blender**

Compact and unique Air / Oxygen blender to maintain perfect FIO₂ from 21% to 100%. Oxygen flow directly can be connected to central pipe line or through flow meter regulator easy to use

- **Air Compressor**

Low Noise Medical grade oil free low pressure Air compressor with bacteria filter, Maximum flow 15L/minute, High pressure relief value protection with constant air flow to the blender

- **Bubble Generator**

CPAP Bubble generator with adjustable CPAP from 3 to 10 cmH₂O, water chamber ease to fill and remove mounted on side of the stand

***Accessories shown in the image may not be part of the standard machine.**