

Infant Incubator-INC18

Product Model : **INC18**



Category: Baby Care Equipment, Medical Equipment, Infant Incubator,

Specification :

- The temp. controlled unit is operated by advanced micro computer technology.
- The three digital display indication provide the Set Temperature, Air Temperature and Skin Temperature.
- Two function control modes provide the clinician more flexibility while optimizing the thermal environment.
- In addition the skin temperature is display is larger than the set / air temperature display.
- The removable heater and control unit make it easy to clean and service.
- Full visibility & accessibility of the baby is provided by the designed of acrylic hood.
- A full complement of alarm provides safety convenience, with both audible and visual indicators :
 - Probe Failure
 - Power Failure
 - High Temperature
 - Low Temperature
- Skin Temperature display can be converted from C to F by push of a button.
- A big LCD is provided for more specific information & visualization of the heater power in different modes.
- The High & Low Temperature alarm setting can be adjusted according to the patient condition.
 - Humidification is provided with the help of two water trays provided at the back of the unit.

- Three side trays provided at the bottom of the unit to keep consumables.
- The unit rests on four big castors with brakes in the front castors.

• **Optional:** Oxygen Delivery System, Electric Slow Suction, CFL Phototherapy

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com