

Infant Incubator-INC18

Product Model : **INC18**



Category: Baby Care Equipment, Medical Equipment, Infant Incubator,

Legacy of Excellence in Healthcare...

Description

- The temp. controlled unit is operated by advanced microcomputer technology.
- The three digital display indications provide the Set Temperature, Air Temperature, and Skin Temperature.
- Two function control modes provide the clinician more flexibility while optimizing the thermal environment.
- In addition, the skin temperature display is larger than the set/air temperature display.
- The removable heater and control unit make it easy to clean and service.
- Full visibility & accessibility of the baby is provided by the design of an acrylic hood.
- A full complement of alarms provides safety convenience, with both audible and visual indicators :
 - Probe Failure
 - Power Failure
 - High Temperature
 - Low Temperature
 - Skin Temperature display can be converted from C to F by the push of a button.
 - A big LCD is provided for more specific information & visualization of the heater power in different modes.
 - The High & Low Temperature alarm settings can be adjusted according to the patient's condition.
 - Humidification is provided with the help of two water trays provided at the back of the unit.
 - Three side trays are provided at the bottom of the unit to keep consumables.
 - The unit rests on four big castors with brakes in the front castors.
- **Optional:**
Oxygen Delivery System, Electric Slow Suction, CFL Phototherapy