

## **Infant Incubator-INC18**

Product Model : **INC18**



**Category: Baby Care Equipment, General Medical Devices, Infant Incubator,**


### **Specification :**

- The temp. controlled unit is operated by advanced micro computer technology.
- The three digital display indication provide the Set Temperature, Air Temperature and Skin Temperature.
- Two function control modes provide the clinician more flexibility while optimizing the thermal environment.
- In addition the skin temperature is display is larger than the set / air temperature display.
- The removable heater and control unit make it easy to clean and service.
- Full visibility & accessibility of the baby is provided by the designed of acrylic hood.
- A full complement of alarm provides safety convenience, with both audible and visual indicators :
  - Probe Failure
  - Power Failure
  - High Temperature
  - Low Temperature
- Skin Temperature display can be converted from C to F by push of a button.
- A big LCD is provided for more specific information & visualization of the heater power in different modes.
- The High & Low Temperature alarm setting can be adjusted according to the patient condition.
  - Humidification is provided with the help of two water trays provided at the back of the unit.

- Three side trays provided at the bottom of the unit to keep consumables.
- The unit rests on four big castors with brakes in the front castors.

• **Optional:** Oxygen Delivery System, Electric Slow Suction, CFL Phototherapy

*Legacy of Excellence in Health Care..*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)