

## **Gynaec/ Obstetric Operation Table, Hydraulic**

Product Model : **OPT107**



**Category: Operation Table, Medical Equipment,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- High quality, reliability and functionality Operating Table is being acknowledged as the ideal solution for Gynecology and Obstetrics.
- GWS-OPT107 is Hydraulic operation table for gynecology and obstetrics.
- Pink color is available for the mattress.

### ***Technical Specifications:-***

- Top Dimension: 1630 mm(L) x 550 mm(W)
- Height Adjustment: 680 mm – 1090 mm
- Seat Plate(down/up): 8 / 15 degree
- Back Plate(down/up): 70 / 10 degree
- Patient Weight Capacity: 200 kgs

Optional:

### **Electric Model (OPT107-E)**

- Battery Backup: 2 hours (optional)
- Voltage: AC 220V $\pm$ 10% 50Hz