


Ophthalmic Operation Table, Hydraulic

Product Model : **OPH13**



Category: Operation Table, Medical Equipment, Ophthalmic Operation Table,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Ophthalmic OT Table has a maximum height of 920 mm and minimum height of 690 mm.
- Ideal for standing table or sitting working position.
- It is an ideal operating table with the advantages of strong applicability, easy to operate.
- Adjusts conveniently and steady performance.
- GWS-OPH13 its ergonomic design ensures both surgeon and patient comfort.
- Patient comfort is assured with cushions which are easy to clean.
- The adjustable cushioned head-rest ensures doctor comfort and resistant during operations.

Technical Specification:-

- Top Dimension: 1900 mm (L) x 600mm (W)
- Width of Table Top 520mm
- Minimum Height: 690 mm
- Maximum Height: 920 mm
- Stroke length: 220 mm
- Patient Weight Capacity: 170 kgs.
- Back rest section: 20? (head side)
- Leg Section: 30?

Optional:-

Electric Model(OPH13-E)

- Power Requirement: 230V AC /50Hz, 110V AC /60 Hz
- Fuses: 2A Slow Blow
- Power in VA: 150VA