

Infant Radiant Warmer with automatic manual & baby(servo)modes

Product Model : **INC19**



Category: Baby Care Equipment, Medical Equipment, Infant Radiant Warmer,

Summary :

GWS brand '**Infant Radiant Warmers**' are overhead heating units that consists of a heat source, a skin-temperature sensor, an automatic/manual baby (servo) control unit, and visual and audible alarms. The device is used to provide thermal support for newborns in the delivery suite, for critically ill infants who require constant nursing intervention, and for infants undergoing treatment that prolongs exposure to a cool environment. With ISO 9001-2015; ISO 13485-2016 certification and the WHO-GMP compliance we can ensure the best quality and medical safety of our product.

Legacy of Excellence in Healthcare...

✉ info@gwsmed.com, gwssurgicals@gmail.com

☎ +91 124 4915548

🌐 www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

