

## Baby Bassinet with utility box

Product Model : **HF-210**



**Category: hospital furniture, General Medical Devices, Baby Crib & Cradles,**

### Specification :

- Frame made up of ERW tubular pipe of 2.5cms. dia. Mounted on four 5cms dia. castors.
- Removable CRIB of poly glass titable to trendelenburg and R-trend positions.
- Laminated bottom shelf with a spacious cabinet for storage.
- Provision for monitor stand and drip rod available.
- Pre-treated and epoxy powder coated.

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com