

Baby Bassinet with utility box

Product Model : **HF-210**



Category: hospital furniture, Medical Equipment, Baby Crib & Cradles,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Frame made up of ERW tubular pipe of 2.5cms. dia. Mounted on four 5cms dia. castors.
- Removable CRIB of poly glass titable to trendelenburg and R-trend positions.
- Laminated bottom shelf with a spacious cabinet for storage.
- Provision for monitor stand and drip rod available.
- Pre-treated and epoxy powder coated.