

## Stretcher on Trolley

Product Model : **HF-133**



**Category: hospital furniture, Medical Equipment, Stretchers,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

- Removable with fixed height.
- Frame is made of CRCA tube.
- Trolley mounted over 150mm castors.
- IV/Saline stand Provision.
- Dimension(approx.): 2100 L x 560 W x 810 H(mm)