

Stretcher on Trolley

Product Model : **HF-133**



Category: hospital furniture, General Medical Devices, Stretchers,

Specification :

- Removable with fixed height.
- Frame is made of CRCA tube.
- Trolley mounted over 150mm castors.
- IV/Saline stand Provision.
- Dimension(approx.): 2100 L x 560 W x 810 H(mm)

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com