

Dressing Trolley (with Bowl & Bucket)

Product Model : **HF-172**



Category: hospital furniture, Medical Equipment,

Summary :

Dressing Trolley is a versatile and durable medical cart designed to facilitate efficient dressing and patient care. Constructed from high-quality stainless steel or other durable materials, it features multiple shelves and compartments for organized storage of medical supplies and instruments. The trolley is equipped with smooth-rolling casters for easy mobility and a sturdy, ergonomic design for user convenience. Ideal for use in hospitals, clinics, and healthcare facilities, the dressing trolley streamlines dressing procedures and supports effective patient care and management.

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Made up of SS tubes mounted on 100mm castors.
- Two SS shelves top with 3 side railing and bottom with 4 side railing.
- Bowl and bucket holding case on one side with SS bowl and bucket.
- Dimension(aprx.): 1010 L x 460 W x 810 H (mm)