

Dressing Trolley (with Bowl & Bucket)

Product Model : **HF-172**




Category: hospital furniture, General Medical Devices,

Specification :

- Made up of SS tubes mounted on 100mm castors.
- Two SS shelves top with 3 side railing and bottom with 4 side railing.
- Bowl and bucket holding case on one side with SS bowl and bucket.
- Dimension(aprx.): 1010 L x 460 W x 810 H (mm)

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com