

Bedside Locker(General)

Product Model : **HF-164**



Category: hospital furniture, General Medical Devices, Bedside Lockers,

Specification :

- Made up of CRCA pipe/sheet with SS top.
- One storage box.
- CRCA tubular leg mounted on PVC stump.
- Pre-treated & epoxy powder coated.
- Dimension (apprx.): 405 L x 405 W x 810 H (mm)

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com