

## Bedside Locker(General)

Product Model : **HF-164**



**Category: hospital furniture, Medical Equipment, Bedside Lockers,**

### Specification :

- Made up of CRCA pipe/sheet with SS top.
- One storage box.
- CRCA tubular leg mounted on PVC stump.
- Pre-treated & epoxy powder coated.
- Dimension (apprx.): 405 L x 405 W x 810 H (mm)

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com