

Obstetrics / Labour Table

Product Model : **HF-136**



Category: hospital furniture, Medical Equipment, Obstetric Labour Tables,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Three section top made up of SS.
- Leg section sliding completely under the main section.
- Base frame made up of CRCA tubes.
- Trendlenburg's position achieved by pneumatic pump.
- Backrest adjustable with ratchet mechanism.
- Lithotomy crutches / Lithotomy rods with straps.
- Pre-treated & epoxy powder coated.
- Dimension(approx.): 860 L x 680 W x 760 H (mm)
- Optional:- Mattress (HF-136.01)