

Double Bowl Stand

Product Model : **HF-189**



Category: hospital furniture, Medical Equipment, Bucket, Basins & stands,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Made with MS tubular frame
- Two bowl fitted in ring of 35cm
- Supplied with two litre stainless steel bowl
- Mounted on 5cm dia. castors
- Pre-treated and epoxy powder coated finish