

Blood Donor Chair

Product Model : **HF-196**



Category: hospital furniture, General Medical Devices, Blood Donor Chairs,

Specification :

- Made up of CRCA pipes/sheet.
- Frame made of MS tubes.
- Three sections adjustable by manually.
- All three sections with cushioned cover.
- Pre-treated & epoxy powder coated.
- Dimension: L 760 x W 800 x H 1100mm

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwsurgicals@gmail.com  +91-124-4915548  www.gwsmed.com