

Blood Donor Chair

Product Model : **HF-196**



Category: hospital furniture, Medical Equipment, Blood Donor Chairs,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Made up of CRCA pipes/sheet.
- Frame made of MS tubes.
- Three sections adjustable by manually.
- All three sections with cushioned cover.
- Pre-treated & epoxy powder coated.
- Dimension: L 760 x W 800 x H 1100mm