

Two Functional Fowler Bed Mechanical Standard (Tubular Panels)

Product Model : **HF-121**



Category: hospital furniture, Medical Equipment, Hospital Beds,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- 4 Section divided mattress area.
- Tubular head & foot bows of unequal height with vertical tube support.
- Saline stand provision on four locations.
- Pre-treated & epoxy powder coated.
- Overall approx size: 2090 L x 910W x 600H (mm).
- Functions by Smooth Crank System: (a)Backrest elevation (b)Kneerest elevation

Optional:-

- Four Section Mattress
- S.S Telescopic IV Rod
- S.S Collapsible Side Railing