

## Two Functional Fowler Bed Mechanical Standard (Tubular Panels)

Product Model : **HF-121**



**Category: hospital furniture, Medical Equipment, Hospital Beds,**

### Specification :

- 4 Section divided mattress area.
- Tubular head & foot bows of unequal height with vertical tube support.
- Saline stand provision on four locations.
- Pre-treated & epoxy powder coated.
- Overall approx size: 2090 L x 910W x 600H (mm).
- Functions by Smooth Crank System: (a)Backrest elevation (b)Kneerest elevation

### Optional:-

- Four Section Mattress
- S.S Telescopic IV Rod

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com