

## HI-VAC, MS Frame

Product Model : **SU-101**



**Category: Suction Machine, Medical Equipment, Electric Suction Units,**

**Legacy of Excellence in Healthcare...**

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)

 +91 124 4915548  [www.gwsmed.com](http://www.gwsmed.com)



**Corp Office:** 1205, 12<sup>th</sup> Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

**Works:** Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

## Description

| Housing                | MS Powder coated on four castors |
|------------------------|----------------------------------|
| Max. Vacuum/LPM        | -700mm Hg, 50 Ltr. / min.        |
| Pump Type              | Double Rotary Vane Pump          |
| Jars                   | 2 x 2 Ltr. PC Jar                |
| Filter                 | Trap Filter                      |
| Tubing                 | 10 mm ID x 2 mtr. (PVC)          |
| Vacuum Gauge           | 10 cm dia. 0-760 mm Hg.          |
| Power Supply           | 220 V AC, 50 / 60 Hz , 380 watt  |
| Noise Level            | 50 dB A $\pm$ 3                  |
| Dimensions             | 49 x 36 x 70 cm                  |
| Net Weight             | 37 Kg                            |
| Overflow Safety        | Mechanical                       |
| Accessories / Optional | 110 V on request                 |