

Electric Cum Foot operated Suction

Product Model : **SU-23**



Vacuum Pump



600ml PC Jar



Oil Free

Category: Suction Machine, Medical Equipment, Electric cum Manual Suction Units,

Legacy of Excellence in Healthcare...

Description

- Housing: Powder Coated anti corrosive MS
- Jars: 600ml Poly-carbonate
- Overflow protection: Mechanical type
- Tubing: Non-collapsible suction tubing.
- Vacuum gauge: 2.0 inch, 0-760mm Hg
- Filter: Autoclavable bacterial filter.
- Pump type: Oil free piston & manual pump
- Capacity: -580mm/ Hg +/- 10 at 12~15 LPM
- Power: 220/ 230V AC / 12V DC
- Backup: 120 minutes