

Perforated Instrument Tray without Cover

Product Model : **HH25**



Category: Stainless Steel Hospital Hollowware, Medical Equipment,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India



Description

- Made with High Quality Stainless Steel.
- Standard Height 2" inches

| Item Code | Sizes |
|-----------|-----------|
| HH25.01 | 18" X 12" |
| HH25.02 | 15" X 12" |
| HH25.03 | 12" X 10" |
| HH25.04 | 12" X 8" |
| HH25.05 | 11" X 7" |
| HH25.06 | 10" X 8" |
| HH25.07 | 9" X 6" |