

3 Section Delivery Table Mechanical

Product Model : **HF301**



Category: Obstetric Labour Tables, hospital furniture, General Medical Devices,

Specification :

- Square tube MS Framework, Pre-treated & epoxy powder coated.
- 3 section cushioned top
- Mounted on rubber shoes.
- Supplied with basin.
- Optional:- Lithotomy Pads / Lithotomy Striups
- Knockdown Construction
- 180L x 60W x 82H cms.
- Available in Stainless Steel (Code):- **HF301.SS**

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com