

3 Section Delivery Table Mechanical

Product Model : **HF301**



Category: Obstetric Labour Tables, hospital furniture, Medical Equipment,

Summary :

The complete phase of child birth can be stressful and overwhelming. Our **3 Section Mechanical Delivery Tables** are designed to ease the trauma of birth by providing comfort and support throughout the stages of delivery. A selection of accessories increases the comfort of both the mother and nursing staff. The key features include a square tube MS Framework which is pre-treated and comes with epoxy powder coating. A 3-section cushioned top to ensure the comfort of the patient under labor. The complete bed is mounted on top of rubber shoes to ensure firm grip and is also supplied with a basin. Lithotomy Pads/Striups are included as optional components. The delivery table is also available in stainless steel (Code: HF301.SS).

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- Square tube MS Framework, Pre-treated & epoxy powder coated.
- 3 section cushioned top
- Mounted on rubber shoes.
- Supplied with basin.
- Optional:- Lithotomy Pads / Lithotomy Striups
- Knockdown Construction
- 180L x 60W x 82H cms.
- Available in Stainless Steel (Code):-
HF301.SS