

Gynaec Examination cum Delivery table, Mechanical

Product Model : **HF302**



Category: Obstetric Labour Tables, hospital furniture, Medical Equipment,

Legacy of Excellence in Healthcare...

 info@gwsmed.com, gwssurgicals@gmail.com

 +91 124 4915548  www.gwsmed.com



Corp Office: 1205, 12th Floor Ocus Quantum Sector 51, Gurugram - 122003 (Hr.) India

Works: Plot No. 38, Sector 7, Sidcul Industrial Area, Haridwar-249403 (Uttarakhand), India

Description

- MS Tubular frame, pre-treated and epoxy powder coated.
- Fixed height
- 3 section cushioned top, framework mounted on rubber shoes.
- Padded knee crutches adjustable.
- Knockdown Construction, Sides of this section are fit with vertical handgrips.
- Optional:- Lithotomy Pads / Lithotomy Striups.
- Supplied with S.S Basin
- Available in Stainless Steel (Code):-
HF302.SS